

GMO Detective



Wetware Guide:

The wetware guide will help you to gather and buy all the biological (wet) materials you need to set up the experiments. This consists of reagents and DNA (controls and primers). These companies usually do not sell to private individuals and the prices go down drastically when buying in bulk, so better to contact us :)

Enzyme mix	1
Primers	2
COX1 : Plant Gene primers	2
GMO Gene primers	3
Making the mastermix	5
1. Enzyme mix	

This is a premade mix consisting of an improved Isothermal Strand Displacing DNA polymerase enzyme along with Nucleotides and a buffer that allow the reaction to happen, all you need to add to this is a primer set and sample.

For simplicity we currently use NEB WarmStart® LAMP Kit (DNA & RNA) E1700:

<https://www.neb.com/products/e1700-warmstart-lamp-kit-dna-rna>

Available worldwide through their Subsidiaries / Distributors.

The small package is enough for 250 reactions which is 30 Experiments.
(a regular LAMP reaction is 25ul, we do 8x10ul reactions per experiment)

2. Primers

The experiment uses 2 sets of QUASR LAMP primers sets specific to are certain region. Each set has 7 Primers :

2 long internal primers

2 short external primers

2 loop primers to make the reaction faster

1 quenching primer to quench the fluorescence when the reaction is negative

Primers can be ordered from you Local Supplier, but are cheaper in bulk so contact us for purchasing together :)

Standard desalination is fine, except for quencher probes which need to be HPLC purified

ldtdna.com

Eurofinsgenomics.com

Thermofisher.com

Sigmaaldrich.com

a. COX1 : Plant Gene primers

COX1 gene common to all plants

Based on: Rapid Detection of *Phytophthora ramorum* and *P. kernoviae* by Two-Minute DNA Extraction

Followed by Isothermal Amplification and Amplicon Detection by Generic Lateral Flow Device J. A.

Tomlinson, M. J. Dickinson, and N. Boonham

doi:10.1094/PHYTO-100-2-0143

Primer	Sequence 5'to3'	Modification
COX F3	tatgggagccgttttgc	
COX B3	aactgctaagrgcattcc	5' FAM
COX FIP	atggattgrcctaaagtctcagggcaggatttcactattgggt	
COX BIP	tgcatttcttagggctttcggatccrgcgtagcatctg	
COX F-Loop	atgtccgaccaaagattttacc	
COX B-Loop	gtatgccacgtcgcattcc	
COX FIPQ	YCAAATCCAT	3' Black Hole Quencher®-1 or Iowa Black® FQ

b. GMO Gene primers

CaMV-35S promoter common to most GMO plants:

Based on: Real-time loop-mediated isothermal amplification for the CaMV-35S promoter as a screening method for genetically modified organisms; Fukuta, S., Mizukami, Y., Ishida, A. et al. Eur Food Res Technol (2004) 218: 496.

<https://doi.org/10.1007/s00217-003-0862-5>

Primer	Sequence 5'to3'	Modification
35S_FIP	aggcatctcaacgatggccttaaaggaaggtggctctc a	
35S_BIP	tgccgacagtggccccaaagtgaagacgtggttggaacg	5' FAM
35S_F3	tgcccagctatctgtcactt	
35S_B3	tcccttacgtcagtgagat	
35S_FLoop	tcctttatcgcaatgatg	
35S_BLoop	agcatcgtggaaaaagaag	
35S_BIPQ	ACTGTCGGCA	3' Black Hole Quencher®-1 Or Iowa Black® FQ

Primer	Primer working stock conc.	Primer conc. in 10µl LAMP rxn	Vol. of primer (1 rxn) µL	Vol. of primer (10 rxn) µL	Vol. of primer (50 rxn) µL	Vol. of primer (100 rxn) µL	Vol. of primer (200 rxn) µL	Vol. of primer (500 rxn) µL
FIP	100µM	1.6µM	0.2	2	10	20	40	100
BIP	100µM	1.6µM	0.2	2	10	20	40	100
F3	100µM	0.2µM	0.025	0.25	1.25	2.5	5	12.5
B3	100µM	0.2µM	0.025	0.25	1.25	2.5	5	12.5
FLP	100µM	0.8µM	0.1	1	5	10	20	50
BLP	100µM	0.8µM	0.1	1	5	10	20	50
Quencher	100µM	2.0µM	0.3	3	15	30	60	150
DNaseFree H2O			2.05	20.5	102.5	205	410	1025
Final volume			3	30	150	300	600	1500

3. Making the mastermix

We want to make enough mastermix for **X** experiments

Each experiment is 8 tubes, 4 for plant gene and 4 for the GMO gene.

We mix $5\mu\text{l} * 4 * \mathbf{X}$ LAMP mix with $3\mu\text{l} * 4 * \mathbf{X}$ Primer mix for each primer set.

It is better to usually make about 10% extra just in case of mistakes and spilling

Enjoy!